

Genius 8080 - Solução total

Código do produto: XX-XXXX



Descrição

O Genius 1053 e o Genius 3070 foram desenvolvidos com tecnologia testada e comprovada, para atender aos requerimentos de gás do LC-MS 8080 da Shimadzu.



Generator Details

- **Tipos de gases:** Nitrogênio
- **Compressor integrado:** Não

Características principais

- **Versátil** - Disponível com ou sem compressor de ar interno
- **Pureza** - Ar zero com concentração de hidrocarboneto (como metano) inferior a 0,1 ppm
- **Duradouro** - Poucos itens consumíveis substituídos garantem baixo custo com manutenção
- **Solução única** - Um único gerador pode atender a várias aplicações
- **Simple instalação** - Necessidade mínima de configuração
- **Econômico** - Mais econômico do que outra fonte de fornecimento de gás
- **Conveniente** - Gás produzido sob demanda, adequando-se a rotina de seu laboratório

Especificações técnicas

Genius 8080 - Solução total

Tipo de gás	Nitrogênio
Tempo de inicialização:	60 mins
Corrente:	0
Nível de ruído:	0

[PEAK Protected]TM

Peak Scientific gas generators define the benchmark in reliability, convenience and performance in laboratories around the world, and come backed with a 12 month warranty. Beyond this period however you can ensure that your investment continues to be **[Protected]** by our comprehensive generator care cover.

Our world-class aftercare support packages deliver a program of scheduled preventative maintenance whilst giving you the reassurance of instant access to worldwide technical support and priority on-site response in the untimely event of a breakdown.

- For ordering parts visit: www.peakscientific.com/ordering/
- For service plans visit: www.peakscientific.com/service/service-plans/

Peak Scientific UK
Tel: +44 (0)141 812 8100
Fax: +44 (0)141 812 8200

Peak Scientific North America
Tel: +1-800-767-6532
Fax: +1-978-608-9503

Peak Scientific China
Tel: +86 21 5079 1190
Fax: +86 21 5079 1191
For a full list of our worldwide office locations, please visit:
Web: www.peakscientific.com
Email: discover@peakscientific.com

Peak Scientific's Quality Management System conforms to: ISO:9001

